
## Uterine Manipulator

### For use in the following procedures:

Diagnostic Laparoscopy, Minilaparotomy, Fertility Examinations, Salpingoplasties  
**Sterile** – (Ethylene Oxide Gas Sterilized) – unless package has been opened or damaged  
**Disposable** – Discard after single use

#### Uterine Manipulator

- |   |                      |
|---|----------------------|
| A. Inflation Valve & Pilot Balloon Assembly | E. Cervical stop     |
| B. Distal Endport                           | F. Handle            |
| C. Cuff                                     | G. Luer Lock Adaptor |
| D. Markings (cm)                            |                      |

#### DESCRIPTION:

The Uterine Manipulator is both a uterine manipulator and a uterine injector for single use. This is a sterile (by ethylene oxide) disposable product made of clear polyvinyl chloride which meets USP recommendations for implant testing. The product is designed with a double lumen, one for inflation of a 10cc intrauterine cuff and the other for injection of fluid through a distal endport. The product is curved to facilitate forward uterine manipulation. The product features an inflation valve and pilot balloon assembly (A), and endport (B), an inflatable cuff (C), centimeter depth markings (D), cervical stop (E), a removable rigid plastic handle (F), and a luer fitting to accommodate a syringe (G). The instrument has a length of 13 inches (33cm) and an outside diameter of 4.5mm (.18 inches).

#### INDICATIONS:

Use of the Uterine Manipulator is indicated in Diagnostic Laparoscopy, Minilaparotomy, Fertility Examinations, and Salpingoplasties procedures where manipulation of the uterus is required. This product also facilitates the sealing of the cervical os while providing a fluid or air injection port.

#### CONTRAINDICATIONS:

- Pregnancy
- Uterine or tubal infection

#### WARNINGS:

- Insure that the handle is secured to the catheter with the flat face of the cervical stop well fixed on the tube at the 6cm mark (see Figure 1). It is possible that the handle may be loosened during shipment or handling or completely off the tube. It is easy to snap back on at the 6cm mark before usage.
- Inflation of cuff before manipulation is essential to avoid uterine trauma.
- Insert the Uterine Manipulator along the proper axis to avoid trauma to the uterus.

#### PRECAUTIONS:

- Lubricate catheter tip before insertion.
- Test inflatable cuff before insertion for possible leakage.
- Check for necessity to dilate cervix before insertion to avoid tearing inflatable cuff.
- When injecting any liquid media adhere to appropriate manufacturer's advice.
- After insertion and inflation insure that inserted cuff (C) is properly inflated by squeezing pilot balloon (A) and checking for tautness. An instrument with a deflated cuff may injure or perforate the uterine wall.

#### ADVERSE REACTIONS:

- Cramping
- Infection
- Perforation of uterine wall

#### DIRECTIONS FOR USE:

1. Check to see that the sterile bag containing the Uterine Manipulator has not been punctured or damaged thus destroying sterility.
2. Be sure handle (F) is securely in place with flat face of cervical stop at 6cm.
3. Using a standard plastic syringe, test the inflatable cuff by inserting 6 to 7cc of air into valve assembly (A).
4. After testing inflation, evacuate all air from the cuff for insertion of the manipulator. With the patient in the lithotomy position, expose the cervical os. Then grasp the anterior lip of the cervix with a single tooth tenaculum.
5. Sound the uterus for depth and direction. (Unless the uterus sounds to less than 6cm, insert the Uterine Manipulator as is with the face of the cervical stop securely attached at the 6cm mark on the tube. The handle is removable for adjustment.)
6. Ascertain if the cervix needs dilation. Dilate to #14 Hank size using approved surgical techniques or use a closed standard uterine dressing forceps to dilate the cervix. The inflatable cuff (C) of the Uterine Manipulator may tear if a wide enough passage is not afforded.
7. Lubricate the distal endport (B) and cuff (C) for easy insertion. Carefully guide the instrument along the natural axis of the cervix to avoid injury.
8. Insert the Uterine Manipulator fully into the uterus until the face of the cervical stop (E) abuts the external cervix. Slowly inflate the cuff with 3 to 10cc of air until you notice the cervical stop being pulled up tightly against the cervix. The inflating cuff "grasps" the lower uterus between itself and the fixed cervical stop. Hold thumb on plunger as you remove syringe from inflation valve to prevent reflux of air back into syringe.
9. To insure that the cuff (C) has not ruptured during procedure, check the tautness of the pilot balloon (A). A soft balloon indicates a ruptured or leaking cuff.
10. Pull gently on instrument to be sure cuff inflation is adequate to prevent easy expulsion of the cuff during injection procedures. Excessive cuff inflation may cause utero-tubal spasm (foreign body reaction) and physiologic closure of patent tubes. Rapid injection of fluid (or gas) or injection of cold fluid may also initiate utero-tubal spasm.
11. The speculum and tenaculum can now be removed since the uterus is secured and sealed. The patient can be placed in the dorsoecumbent position with the proximal end of the Uterine Manipulator between the legs.
12. Posterior positioned uteri (see Figure 2) can be rotated to the anterior position by placing the Uterine Manipulator into the uterus posteriorly and rotating it anteriorly after cuff inflation.
13. Before manipulating or injecting, check that the cuff (C) is inflated by observing pilot balloon inflation. The cervical seal will be lost should cuff break. A deflated cuff increases the danger of uterine wall perforation.
14. To remove the Uterine Manipulator insert syringe deeply in the air valve (A) and deflate cuff. Carefully remove the instrument. Check before discarding instrument to be sure it is intact and no parts remain in the vaginal canal.


Figure 2

#### PRODUCT INFORMATION:

Product Number: TM1151  
 OD Size: 4.5mm  
 Length: 33cm (13")  
 Packaging: 12 sterile units per box

#### MANUFACTURED BY:

**Catheter Research, Inc.**  
 dba Thomas Medical, Inc.  
 5610 West 82nd Street  
 Indianapolis, IN 46278 USA  
 Telephone: (800) 556-0349  
 (317) 872-0074  
 Facsimile: (317) 872-0169

Emergo Europe  
 Molenstraat 15  
 2513 BH The Hague  
 The Netherlands  
 Telephone: 31(0) 70 345-8570  
 Facsimile: 31(0) 70 346-7299

Comments regarding this device can be directed to Attn: Quality Assurance Dept., CRI, 5610 West 82nd Street, Indianapolis, IN 46278 USA, Fax: (317) 872-0169

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