Uterine Injector-4.0mm

For use in the following procedures:

Rubin's Test, Hysterosalpingograms, and Salpingoplasties

Sterile – (Ethylene Oxide Gas Sterilized) – unless package has been opened or damaged

Disposable – Discard after single use


Figure 1

Uterine Injector-4.0mm

- A. Inflation Valve & Pilot Balloon Assembly
- B. Distal Endport
- C. Cuff
- D. Markings (cm)
- E. Cervical stop
- F. Stopcock

DESCRIPTION:

The **Uterine Injector-4.0mm** is a single use, sterile/disposable, clear polyvinyl chloride uterine injector which meets USP recommendations for implant testing. This product is designed with a double lumen, one for inflation of a 5cc intrauterine cuff and the other for injection of fluid through a distal endport. The product is slightly curved to facilitate easy entrance into the uterus through the cervical canal. The product features (see Figure 1) an inflation valve and pilot balloon assembly (A), an endport (B), an inflatable cuff (C), centimeter depth markings (D), a removable cervical stop (E), and a stopcock to accommodate a syringe (F). The device has a length of 23cm and an outside diameter of 4.0mm.

INDICATIONS:

The **Uterine Injector-4.0mm** is indicated for procedures such as Rubin's Test, hysterosalpingography, salpingoplasties, or hydrotubation. It is designed for "in-office" use. The **Uterine Injector-4.0mm** can be used without cervical dilation or anesthesia if the uterus can be sounded easily with a standard uterine sound. This device provides for efficient sealing of the cervix when injections of gas or liquid are required.

CONTRAINDICATIONS:

- Pregnancy
- Uterine or tubal infection

WARNINGS:

- Insure that the snap-on cervical stop (E) is secured to the catheter with the flat face of the cervical stop well fixed on the tube at the 6cm mark (see Figure 1). It is possible that the stop may become loosened from shipment or improper handling. It is easy to snap back on at the 6cm mark *before* usage.
- To avoid uterine trauma, insert the **Uterine Injector-2.0mm** along the natural axis of the cervical canal.

PRECAUTIONS:

- Inflate uterine cuff (C) to check for leakage *before* insertion.
- Lubricate distal end (cuff and catheter) for easy insertion.
- Observe manufacturer's directions and precautions pertaining to any liquid media or gas that is being injected into the patient.

ADVERSE REACTIONS:

- Cramping
- Infection
- Perforation of uterine wall

INSTRUCTIONS FOR USE:

1. Check to see that the sterile pouch containing the **Uterine Injector-4.0mm** has not been punctured or damaged thus compromising sterility.
2. Inspect carefully to ensure that the cervical stop (E) is securely in place with the face at the 6cm mark. If the stop is not fixed and slides easily, snap it securely in place before using.
3. Using the plastic syringe provided, inflate the cuff (C) with 5cc of air to check for leakage before inserting into the patient. Air is injected through the valve assembly (A).
4. After checking inflatability, evacuate all the air.
5. Using a speculum, expose cervix and grasp the anterior lip with a single tooth tenaculum.
6. Sound the uterus for depth and direction. If the sound does not slide easily or drags, do not force the injector through. The cuff may tear if the cervix opening is too narrow.
7. Lubricate the distal endport (B) and cuff (C) for easy insertion. Carefully guide the instrument along the natural axis of the cervix to avoid injury.
8. Insert the **Uterine Injector-4.0mm** fully into the uterus until the face of the cervical stop (E) abuts the external cervix. Slowly inflate the cuff with 3 to 5cc of air until the cervical stop is pulled up tightly against the cervix. The inflating cuff "grasps" the lower uterus between itself and the fixed cervical stop. Hold thumb on plunger as syringe is removed from the inflation valve to prevent reflux of air back into the syringe.
9. To insure that the cuff (C) has not ruptured during procedure, check the tautness of the pilot balloon (A). A soft balloon indicates a ruptured or leaking cuff.
10. Pull on the instrument gently to see if the inflation of the cuff is adequate to prevent expulsion of the inflated cuff with normal intrauterine injection pressure. The less the cuff is inflated, the easier it can be expelled with excessive intrauterine injection pressure. (If, however, the cuff is inflated beyond the amount of air needed to gently grasp the uterus, there will be more discomfort for the unanesthetized patient. The tense inflated cuff in the smaller uterus may also set up a "foreign body" like reaction causing uterine spasm that may give physiological blockage of the tubes at the cornua).
11. The uterus is now "grasped" and the os sealed. The speculum and tenaculum can now be removed without disturbing the seal or losing the device. The patient can be placed in the dorsorecumbent position with easy access to the uterine cavity maintained. The plastic tube is pliable enough to be bent in the desired direction.
12. Injection media may be applied with a standard syringe or extenders as needed. Easy access to the uterus during salpingoplasties may be achieved by attaching an extender to the stopcock (F).
13. To remove the **Uterine Injector-4.0mm**, insert the syringe deeply into the inflation valve (A) and deflate the cuff. Carefully remove the instrument, checking carefully to be sure the device is intact and no parts remain in the vaginal canal.

PRODUCT INFORMATION:

Product Number: TMI1152
OD Size: 4.0mm
Length: 23cm (9")
Packaging: 12 sterile units per box

MANUFACTURED BY:

Catheter Research, Inc.
dba Thomas Medical, Inc.
5610 West 82nd Street
Indianapolis, IN 46278
Telephone: (800) 556-0349
(317) 872-0074
Facsimile: (317) 872-0169

Emergo Europe
Molenstraat 15
2513 BH The Hague
The Netherlands
Telephone: 31(0).70.345.8570
Facsimile: 31(0).70.346.7299

Comments regarding this device can be directed to Attn: Quality Assurance Dept., CRI, 5610 West 82nd Street, Indianapolis, IN 46278 USA, Fax: (317) 872-0169

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