**INSTRUCTIONS FOR USE**  
**ENDOMETRIAL SAMPLER (Endometrial Suction Curette)**  
 Not made with natural rubber latex.

**For use in the following procedures:**

Historical Biopsy of the uterine mucosal lining or sample extraction of uterine menstrual content.

**Figure 1**

**Endometrial Sampler**

- A. Curette Opening
- B. Plastic Sheath
- C. Piston Rod

**PRECAUTIONS:**

Contents Sterile unless package has been opened or damaged.

Single Patient Use Only – Discard after using.

**DESCRIPTION:**

The Endometrial Sampler is a single use, sterile, disposable, suction curette for obtaining a histologic biopsy of the uterine mucosal lining or sample extraction of uterine menstrual content for microscopic examination or culturing.

The device consists of a clear flexible, plastic sheath that is 25cm in length with a 3.0mm OD (outer diameter) and a 2.4mm ID (inner diameter). The sheath is marked with graduation marks at 6cm, 7cm, 8cm, 9cm, 10cm and 11cm distance from the distal tip of the sheath to indicate the depth of insertion of the sheath into the uterus during use.

The extreme distal tip of the sheath is tapered. At the side of the distal tip of the sheath is a curette opening which is 2.0mm in diameter, has sharp edges and leads to the lumen of the sheath.

A piston rod can be moved forward and backward with the full length of the lumen of the sheath. This is accomplished by digital manipulation of the knob affixed to the proximal end of the piston rod which extends beyond the proximal end of the sheath. The piston is prevented from being totally pulled from within the sheath by means of an indentation in the sheath located at its proximal end.

Rapid movement of the piston within the sheath from its fully inserted position to its maximum retracted position creates a negative pressure (suction) within the lumen of the sheath. This negative pressure draws the mucosal tissue through the curette opening and into the lumen of the sheath as the curette scrapes against that endometrial wall while it is maneuvered within the uterine cavity.

**INDICATIONS:**

Histologic biopsy of the glandular epithelium and superficial chronic layers of the uterine endometrial wall or simple extraction of uterine menstrual content for:

- Detection of endometrial carcinoma and precancerous conditions that could make estrogen therapy advisable
- Determination of endometrial tissue response to hormonal influences by visual observation of the histologic sample such as:
  - Endometrial dating
  - Periodic scrutiny of the effects of hormonal therapy on the uterine endometrium, particularly in patients receiving estrogen replacement therapy for menopausal symptoms and retardation of osteoporosis
  - Luteal insufficiency (visually determined histologically from endometrial tissue obtained only during the secretory phase of the menstrual cycle)
  - Functional metrorrhagia (contraception, pre-menopause)
- Diagnosis of pathology associated with infertility.

**CONTRAINDICATIONS:**

- Should not be used in patients who are pregnant or suspected of being pregnant.
- Should not be used in the presence of chronic cervical infection.
- Should not be used in the presence of or after recent pelvic inflammatory disease.

**WARNINGS:**

In patients with amenorrhea, obtaining an endometrial biopsy with the Endometrial Sampler should be performed only after confirmation of the absence of detectable circulating HCG levels.

**PRECAUTIONS:**

- By virtue of its flexibility and small caliber, with an Endometrial Sampler a histologic endometrial biopsy can be performed with minimal danger of uterine wall perforation. However, care should be taken prior to its insertion to ascertain the depth of the uterus and any uterine retroflexion or antifixion present in the individual patients. Best clinical judgement should be used to allow for any variations from anatomical norms during insertion and in no case should the instrument be forced against digitally felt resistance.
- If passage of the Endometrial Sampler through the cervical canal is not found to be easily achievable, the instrument should not be forced and the possibility of the presence of pathologic cervical stenosis considered. With elderly patients with whom an extremely dry cervix may cause resistance to the instrument's passage, very slight lubrication of the sheath with a water-soluble gel may facilitate its insertion.

**ADVERSE REACTIONS:**

- Perforation of uterine wall
- Cramping
- Uterine spasm

**INSTRUCTIONS FOR USE:**

1. Prepare the vagina and cervix using currently accepted antiseptic techniques as you would for any intrauterine procedure.
2. With a vaginal speculum in place and open to expose the uterine cervix, gently probe the uterus for depth and direction with a uterine sound. Then grasp the cervix with very fine forceps or a tenaculum. If the uterus is anteverted, this grasp should be applied to the anterior lip of the cervix. If the uterus is retroverted, the grasp should be applied to the posterior lip of the cervix. If the uterus is normally positioned and the cervix is stable, it is possible that the use of forceps or a tenaculum may not be necessary.
3. Now withdraw the speculum approximately 2cm and apply gentle traction to the forceps or tenaculum to straighten any cervical curvature and to stabilize the cervix.
4. With its piston positioned at the extreme distal end of the sheath, the Endometrial Sampler is then inserted into and gently passed through the cervical canal and into the cavity of the uterus to a depth that corresponds to that which was determined by uterine depth finding. Slight friction can usually be felt with the fingers as the Endometrial Sampler passes through the internal cervical os. In no event should force be used to obtain passage if resistance is encountered. Slight lubrication of the sheath with a water-soluble gel may facilitate its passage in patients with an extremely dry or narrow cervical canal. If, as is usually not the case, negotiation of the cervical canal proves painful, brief topical use of an anaesthetic agent may help facilitate introduction of the Endometrial Sampler.
5. When the sheath is in position within the uterine cavity, discontinue any traction you have applied with a tenaculum or forceps. Then, while holding the sheath in position with one hand, with the other hand rapidly pull the piston firmly and without interruption (with one swift steady motion) toward the proximal end of the sheath as far as it will go. Slow, interrupted or partial withdrawal of the piston will not produce the amount of negative pressure needed.
6. A slight crimp in the sheath will prevent total withdrawal of the piston from within the sheath. This action creates a negative pressure within the sheath that draws the tissue into the curette opening at the distal end of the sheath and captures the separated tissue within its lumen.
7. Immediately after pulling the piston the full permissible distance the sheath should be rotated between the fingers and moved from side to side simultaneously with 3 or 4 gentle back and forth movements of the instrument within the uterine cavity. It should then be gently withdrawn from the uterus. Examination of the sheath should then reveal a specimen of the uterine mucosa of histologic quality clearly visible within the sheath. Bleeding, if it occurs at all, is usually minimal.
8. For histologic assessment of the specimen, the lip of the sheath should now be sectioned just proximal to its distal curette opening. The piston should then be advanced within the sheath to express the specimen from the sectioned sheath into an appropriate transport medium.


**PRODUCT INFORMATION:**

Product Number: TM11175  
 OD Size: 3.0mm  
 Length: 25cm  
 Packaging: 25 sterile units per box

**MANUFACTURED BY:**  
 Catheter Research, Inc.  
 dba Thomas Medical, Inc.  
 6102 Victory Way  
 Indianapolis, IN 46278  
 Telephone: (800) 556-0349  
 (317) 872-0074  
 Facsimile: (317) 872-0169

**AUTHORIZED REPRESENTATIVE:**  
 Emergo Europe  
 Molenstraat 15  
 2513 BH The Hague  
 The Netherlands  
 Telephone: 31(0).70.345.8570  
 Facsimile: 31(0).70.346.7299

Comments regarding this device can be directed to Attn: Quality Assurance Dept., CRI, 6102 Victory Way, Indianapolis, IN 46278 USA, Fax: (317) 872-0169

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