Sonolject Catheter
for Sonohysterography


Figure 1

Sonolject Catheter™

- A. Luer
- B. Plastic Sheath
- C. Shaft

PRODUCT NO. TMI1170

Device Description:

The Sonolject Catheter consists of a non-balloon catheter and insertion/placement sheath.

Intended Use/Indications:

For administering saline during Sonohysterography procedures to detect uterine pathology such as polyps, fibroids, adhesions or endometrial thickening.

Contraindications:

Suspected infection, suspected pregnancy, profuse bleeding or sexually transmitted disease.

How Supplied:

Sterile in unopened, undamaged package.
Single use only.
Box of 10.

Instructions for Use:

1. Remove catheter from pouch.
2. Attach a saline filled syringe (not supplied) to the luer connector. Inject saline through catheter to remove air. Advance the insertion/placement sheath so that the distal end of the catheter protrudes slightly from the distal end of the sheath.
3. View the cervix and advance the sheath and catheter so that the tip of the catheter enters the cervix.
4. Advance the catheter through the cervical canal and into the uterus.
5. Insert the transvaginal transducer in the proper position adjacent to the catheter. After activation of ultrasound equipment, inject saline into the uterus as required for the procedure. Typically, 3 to 20cc of saline is sufficient for proper diagnosis.
7. Complete the examination and withdraw the catheter.

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Catheter Research, Inc.
dba Thomas Medical, Inc.