K-UMI MANIPULATOR

FEATURES:

- | | |
|--------------------|-----------------------|
| A. Inflation Valve | F. Sliding Handle |
| B. Endport | G. Luer Lock Adaptor |
| C. Balloon | H. Flexible Extension |
| D. Depth Markings | I. Pinch Clamp |
| E. Cervical Stop | J. Pilot Balloon |

PRODUCT INFORMATION:

Product Number: TMI6003
OD Size: 5mm
Length: 33cm


Figure 1

DESCRIPTION:

The **K-UMI Manipulator™** is both a uterine manipulator and a uterine injector for single use. This is a sterile (by ethylene oxide) disposable product made of clear polyvinyl chloride which meets USP recommendations for implant testing. The product is designed with a double lumen, one for inflation of a 10cc intrauterine balloon and the other for injection of fluid through a distal endport. The product is curved to facilitate forward uterine manipulation. The product features an inflation valve and pilot balloon assembly (A), an endport (B), an inflatable balloon (C), centimeter depth markings (D), cervical stop (E), a movable rigid plastic handle (F) with a pinch clamp (I), and a luer fitting to accommodate a syringe (G). The instrument has a length of 33cm (13 inches) and an outside diameter of 5mm.

The K-UMI Manipulator features a sliding handle that can be easily positioned along the length of the catheter shaft. This allows positioning the cervical stop at the desired measured depth. The balloon is inflated in the uterus and retracted against the internal os. The cervical stop is pressed firmly against the cervix by advancing the handle while holding the catheter shaft. When the handle is locked, the K-UMI maintains a secure, trauma-free hold on the uterus. The instrument has a length of 14.5 inches (37cm) and an outside diameter of 5mm (15 Fr).

INDICATIONS:

- Manipulation of the uterus during Laparoscopy, Minilaparotomy, Fertility Examination, and Salpingoplasty.
- Sealing of the cervical os to inject fluid into the uterus.

CONTRAINDICATIONS:

- Pregnancy
 - Uterine or tubal infection

WARNINGS:

- Always position handle at correct sounded depth, using the cm markings on the catheter.
- Inflate the balloon before manipulation to avoid uterine trauma.
- Insert the Uterine Manipulator along the proper axis to avoid uterine trauma.

PRECAUTIONS:

- Test inflatable balloon for possible leakage prior to product use.
- Lubricate catheter tip before insertion.
- Check for necessity to dilate cervix to ease product insertion and to avoid tearing balloon.
- When injecting any liquid media adhere to appropriate manufacturer's advice.

ADVERSE REACTIONS:

- Cramping
- Infection
- Uterine trauma or perforation
- Uterine spasm

DIRECTIONS FOR USE (REF: Figure 1):

- Assure that the package containing the product has not been punctured or damaged.
- Using a standard syringe, test the inflatable Balloon (C) by injecting 6cc of air into Inflation Valve (A). Deflate balloon after test.
- Sound the uterus for depth and direction.
- Squeeze the pinch clamp (I) to position proximal end of the Sliding Handle (F) on catheter shaft at the same cm Depth Marking (D) as sounded depth. If the fundus is sounded at 7 cm, set the proximal end of the handle, indicated by the arrow, at the 7 cm Depth Marking (D). NOTE: A Depth Marking (D) of 7 cm, for example, is equivalent to a distance of 6 cm from Endport (B) to the Cervical Stop (E). This 1 cm difference provides a margin of safety in depth of insertion
- Dilate cervix to #14 Hank size, if necessary.
- Lubricate the Endport (B) and Balloon (C) for ease of insertion. Carefully guide the K-UMI Manipulator along the natural axis of the cervix to avoid injury.
- Advance the instrument into the uterus until the Cervical Stop (E) presses on the cervix. With a syringe filled with 3 to 10cc of air, inflate the balloon through the inflation valve (A). Maintain light pressure to ensure contact between the cervical stop and the cervix while inflating the balloon. Hold thumb on plunger as you remove syringe from Inflation Valve (A) to prevent reflux of air back into syringe.
- Withdraw the catheter slowly until the balloon fits snug against internal OS. Squeeze pinch clamp to advance the cervical stop (E) on the handle (F) against the cervix. Release pinch clamp to lock handle in place. Note the change in Depth Markings (D). The K-UMI Manipulator is now locked and firmly grasping the uterus.
- To insure that the Balloon (C) has not ruptured during procedure, check the tautness of the Pilot Balloon (J). A soft balloon indicates a ruptured or leaking balloon.
- Apply gentle traction on the instrument to be sure it is properly placed.
- Before manipulating or injecting, check that the Balloon (C) is inflated by observing Pilot Balloon (J) inflation. The cervical seal will be lost should balloon break. A deflated balloon increases the danger of uterine trauma.
- Inject through the Luer Lock Adaptor (G) by attaching a syringe. The Flexible Extension (H) provides strain relief for syringe movement.
- To remove the K-UMI Manipulator, insert syringe deeply in the Inflation Valve (A) and deflate Balloon (C). Carefully withdraw the instrument. Check before discarding instrument to be sure it is intact and no parts remain in the patient.

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Catheter Research, Inc.
Db Thomas Medical, Inc.