Sonolject Catheter
for Sonohysterography


Figure 1

Sonolject Catheter™

- A. Luer
- B. Plastic Sheath
- C. Shaft

PRODUCT NO. TMI1170

Device Description:

The Sonolject Catheter consists of a non-balloon catheter and insertion/placement sheath.

Intended Use/Indications:

For administering saline during Sonohysterography procedures to detect uterine pathology such as polyps, fibroids, adhesions or endometrial thickening.

Contraindications:

Suspected infection, suspected pregnancy, profuse bleeding or sexually transmitted disease.

How Supplied:

Sterile in unopened, undamaged package.
Single use only.
Box of 10.

Instructions for Use:

1. Remove catheter from pouch.
2. Attach a saline filled syringe (not supplied) to the luer connector. Inject saline through catheter to remove air. Advance the insertion/placement sheath so that the distal end of the catheter protrudes slightly from the distal end of the sheath.
3. View the cervix and advance the sheath and catheter so that the tip of the catheter enters the cervix.
4. Advance the catheter through the cervical canal and into the uterus.
5. Insert the transvaginal transducer in the proper position adjacent to the catheter. After activation of ultrasound equipment, inject saline into the uterus as required for the procedure. Typically, 3 to 20cc of saline is sufficient for proper diagnosis.
7. Complete the examination and withdraw the catheter.

PRODUCT INFORMATION:

Product Number: TMI1170
OD Size: 5Fr.
Length: 35.1cm
Packaging: 10 sterile units per box

MANUFACTURED BY:
CRI
dba Thomas Medical, Inc.
6102 Victory Way
Indianapolis, IN 46278
Telephone: (800) 556-0349
(317) 872-0074
Facsimile: (317) 872-0169

AUTHORIZED REPRESENTATIVE:
CRI Medical Devices Limited
Westford Enterprise Centre
Strandfield Business Park
Kerlogue, Rosslare Road
Westford, Ireland

Comments regarding this device can be directed to Attn: Quality Assurance Dept., CRI, 6102 Victory Way, Indianapolis, IN 46278 USA, Fax: (317) 872-0169

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