



LED DISPOSABLE SPECULUM

SEE PRODUCT LABEL FOR:

STERILE EO



LED Disposable Speculum - English
LEDSPECU-IF1

MEDGYN PRODUCTS, INC.

100 W. Industrial Rd. Addison, IL 60101 USA

e: info@medgyn.com | www.medgyn.com

t: 630.627.4105

f: 630.627.0127

toll-free: 800.451.9667

DEVICE DESCRIPTION:

MedGyn LED Speculums are made with clear polystyrene for easy, single-hand operation and insertion. The Disposable LED Speculum offers bright, white light to improve visualization of the exam area. The LED light is cold, so it will not increase patient discomfort.

INDICATIONS FOR USE:

Intended for use with vaginal examinations and/or the taking of cervical smears.

CAUTION:

Federal law restricts this device to sale by the order of or for use by a physician.

WARNINGS:

- Do not use if package is open or damaged, or if the seals are not intact.
- This speculum is for single use only. Do not reuse.
- This speculum is disposable and must be discarded after use.
- The device is sterilized by ethylene oxide.

INSTRUCTIONS FOR USE:

1. Open the speculum pouch and, if required, lubricate the speculum with a water soluble lubricant.
2. Ensure the locks are disengaged.
3. With bills in the closed position, part the labia and introduce the speculum gently into the vagina. Once the speculum is in the required position apply gentle pressure to the angle and height adjusters. The speculum will automatically lock into position.
4. To activate the LED light, use the switch on the bottom of the speculum handle.
5. Carry out the procedure as per recommended techniques.
6. On completion of the procedure, withdraw the speculum a short distance with the bills in the open position to avoid entrapment of the cervix. Apply gentle pressure to the locks to disengage. Once the speculum has been withdrawn beyond the cervix then allow the natural pressure of the vaginal walls close the bills and withdraw the speculum slowly.
7. Dispose of the speculum.