## Uterine Injector-2.0mm

For use in the following procedures:

Hysterosalpingograms, Salpingoplasties, and Hydrotubation

**Sterile** – (Ethylene Oxide Gas Sterilized) – unless package has been opened or damaged

**Disposable** – Discard after single use


Figure 1

### Uterine Injector-2.0mm

- A. Inflation Valve
- B. Distal Endpoint
- C. Cuff
- D. Cervical stop
- E. Stopcock
- F. Syringe (not pictured)

### DESCRIPTION:

The **Uterine Injector-2.0mm** is a single use, sterile/disposable, clear polyvinyl chloride uterine injector which meets USP recommendations for implant testing. This product is designed with a double lumen, one for inflation of a 2cc intrauterine cuff and the other for injection of fluid through a distal endpoint. The product is slightly curved to facilitate easy entrance into the uterus through the cervical canal. The product features (see Figure 1) an inflation valve (A), an endpoint (B), an inflatable cuff (C), a removable cervical stop (D), and a stopcock (E) to accommodate a syringe (F). The device has a length of 23cm (9") and an outside diameter of 2.0mm.

### INDICATIONS:

The **Uterine Injector-2.0mm** is indicated for in-office or hospital use when efficient sealing of the uterine cervix is required for the injection of liquid or gas such as hysterosalpingography, salpingoplasties, hydrotubation, or Rubin's Test. The **Uterine Injector-2.0mm** can be used without cervical dilation or anesthesia if the uterus can be sounded easily with a standard uterine sound.

### CONTRAINDICATIONS:

- Pregnancy
- Uterine or tubal infection

### WARNINGS:

- Insure that the snap-on cervical stop (D) is secured to the catheter with the flat face of the cervical stop well fixed on the tube (see Figure 1). It is possible that the stop may become loosened from shipment or improper handling. It is easy to snap back on *before* usage.
- To avoid uterine trauma, insert the **Uterine Injector-2.0mm** along the natural axis of the cervical canal.

### PRECAUTIONS:

- Inflate uterine cuff (C) to check for leakage *before* insertion.
- Lubricate distal end (cuff and catheter) for easy insertion.
- Observe manufacturer's directions and precautions pertaining to any liquid media or gas that is being injected into the patient.

### ADVERSE REACTIONS:

- Cramping
- Infection
- Perforation of uterine wall

### DIRECTIONS FOR USE:

1. Check to see that the sterile pouch containing the **Uterine Injector-2.0mm** has not been punctured or damaged thus compromising sterility.
2. Inspect carefully to ensure that the cervical stop (D) is securely in place. If the stop is not fixed and slides easily, snap it securely in place before using.
3. Using the plastic syringe provided, inflate the cuff (C) with 2cc of air to check for leakage before inserting into the patient. Air is injected through the inflation valve (A).
4. After checking inflatability, evacuate all the air.
5. Using a speculum, expose cervix and grasp the anterior lip with a single tooth tenaculum.
6. Sound the uterus for depth and direction. If the sound does not slide easily or drags, do not force the injector through. The cuff may tear if the cervix opening is too narrow.
7. Lubricate the distal endpoint (B) and cuff (C) for easy insertion. Carefully guide the instrument along the natural axis of the cervix to avoid injury.
8. Insert the **Uterine Injector-2.0mm** fully into the uterus until the face of the cervical stop (D) abuts the external cervix. Slowly inflate the cuff with 1.5 to 2cc of air until you notice the cervical stop being pulled up tightly against the cervix. The inflating cuff "grasps" the lower uterus between itself and the fixed cervical stop. The inflating cuff is removed from the uterus by the inflation valve to prevent reflux of air back into the syringe.
9. Pull on the instrument gently to see if the inflation of the cuff is adequate to prevent expulsion of the inflated cuff with normal intrauterine injection pressure. The less the cuff is inflated, the easier it can be expelled with excessive intrauterine injection pressure. (If, however, the cuff is inflated beyond the amount of air needed to gently grasp the uterus, there will be more discomfort for the unanesthetized patient. The tense inflated cuff in the smaller uterus may also set up a "foreign body" like reaction causing uterine spasm that may give physiological blockage of the tubes at the cornua).
10. The uterus is now grasped and the os sealed. The speculum and tenaculum can now be removed without disturbing the seal or losing the device. The patient can be placed in the dorsorecumbent position with easy access to the uterine cavity maintained. The plastic tube is pliable enough to be bent in the desired direction.
11. Injection media may be applied with a standard syringe and plastic stopcock or extenders as needed. Easy access to the uterus during salpingoplasties may be achieved by attaching an extender to the stopcock (E).
12. To remove the **Uterine Injector-2.0mm**, insert the syringe deeply into the inflation valve (A) and deflate the cuff. Carefully remove the instrument, checking carefully to be sure the device is intact and no parts remain in the vaginal canal.

### PRODUCT INFORMATION:

Product Number: TMI1153  
OD Size: 2.0mm  
Length: 23cm (9")  
Packaging: 12 sterile units per box

### MANUFACTURED BY:

**CRI**  
**dba Thomas Medical, Inc.**  
6102 Victory Way  
Indianapolis, IN 46278  
Telephone: (800) 556-0349  
(317) 872-0074  
Facsimile: (317) 872-0169

**CRI Medical Devices Limited**  
Wexford Enterprise Centre  
Strandfield Business Park  
Kerlogue, Rosslare Road  
Wexford, Ireland

Comments regarding this device can be directed to Attn: Quality Assurance Dept., CRI, 6102 Victory Way, Indianapolis, IN 46278 USA, Fax: (317) 872-0169

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