

FEATURES

25 different sizes ranging from 5mmx5mm to 25mmx10mm
Loop, ball and square electrodes

TECHNICAL DATA

Shaft Length 13cm (5")
Sterilization Method ETO
Box Qty. 5 ea.

ORDERING INFORMATION

LE-29-505 5mm x 5mm, Square Electrode	LE-26-150 15mm x 10mm Loop Electrode
LE-01-104 10mm x 4mm, Square Electrode	LE-08-153 15mm x 13mm Loop Electrode
LE-06-108 10mm x 8mm, Square Electrode	LE-08-151 15mm x 15mm Loop Electrode
LE-06-101 10mm x 10mm, Square Electrode	LE-09-208 20mm x 8mm Loop Electrode
LE-03-505 5mm x 5mm Loop Electrode	LE-24-201 20mm x 10mm Loop Electrode
LE-04-105 10mm x 5mm Loop Electrode	LE-28-212 20mm x 12 mm Loop Electrode
LE-22-107 10mm x 7mm Loop Electrode	LE-10-203 20mm x 13mm Loop Electrode
LE-28-108 10mm x 8mm Loop Electrode	LE-11-205 20mm x 15mm Loop Electrode
LE-05-101 10mm x 10mm Loop Electrode	LE-25-202 20mm x 20mm Loop Electrode
LE-05-103 10mm x 13mm Loop Electrode	LE-12-251 25mm x 10mm Loop Electrode
LE-23-155 15mm x 5mm Loop Electrode	LE-13-003 3mm Diameter Ball Electrode
LE-27-157 15mm x 7mm Loop Electrode	LE-14-005 5mm Diameter Ball Electrode
LE-07-158 15mm x 8mm Loop Electrode	

Download IFU



Thomas Medical is a distributor for all electrosurgical units and accessories.



PREMIER ELECTROSURGICAL UNIT

- Model 2001e
- Output Modes Cut, Coagulation, Blend, Fulguration and Bipolar
- Max Watt Output 120
- Controls Hand or Foot
- Display Performance Dual Digital Displays

ORDERING INFORMATION

- 9006501 Electrosurgical Generator
- Includes: Generator (100-120v), Power Cord, Single Cell Disposable Patient Pad, Disposable Hand Control Pencil, Operating Manual.




SINGLE & DUAL CELL DISPERSIVE PATIENT PAD

- 9006540 Single Cell Pad, 5/pk
- 9006550 Dual Cell Pad, 5/pk



DISPOSABLE PENCILS

- Sold in Box of 20
- LQ-05-020 Hand Switching Pencil and cord with 3 Prongs
- LQ-09-020 Foot Control Handle and Cord



COATED INSTRUMENTS

Please download our catalog to view a complete list of coated instruments. [Click Here to Download](#)

- a. LL-34-450 Large Pederson Coated Speculum with Smoke Tube (4.5"x1")
- b. LL-38-425 Large Collin Speculum with Smoke Tube (4.25"x1.5")
- c. LL-44-100 Schroeder Tenacula with Exposed Tips (25cm)