

**NEW DESIGN**


Every sampler is vacuum tested for reliable performance.



**FEATURES**

- Depth markings guide insertion depth and aid in catheter placement
- Atraumatic tip reduces discomfort upon insertion
- Firm and flexible cannula prevents kinking
- New handle design creates a better grip during retrieval of the sample

**Endometrial sampler is now non-DEHP & PVC Free**

 Not made with natural rubber latex.

**TECHNICAL DATA**

Depth Mark Range . . . . .	5-11 cm	100% Vacuum Tested
Length . . . . .	24.5cm	French Size . . . . . 9.5Fr
Sterilization Method . . . . .	ETO	Box Qty. . . . . 25 ea.

**ORDERING INFORMATION**

TMI1175 . . . . . Endometrial Sampler

Download IFU

