

Everything is in one sterile field to assist in an efficient and effective procedure



### TRAY CONTENTS

- |  |       |   |       |
|--|-------|---|-------|
| 1. Clear Tray 16"x6"x2" . . . . .          | 1 ea. | 9. Sanitary Napkin . . . . .              | 1 ea. |
| 2. Clear Tray 5"x4"x2" . . . . .           | 1 ea. | 10. 3"x3" Gauze Pads . . . . .            | 5 ea. |
| 3. Povidone-Iodine Ointment . . . . .      | 1 ea. | 11. 18-gauge Safety Needle . . . . .      | 1 ea. |
| 4. Povidone-Iodine Solution . . . . .      | 2 ea. | 12. 8" Swabs . . . . .                    | 3 ea. |
| 5. Lubricating Jelly . . . . .             | 1 ea. | 13. 20cc Syringe . . . . .                | 1 ea. |
| 6. 30"x30" Over-wrap (not shown) . . . . . | 1 ea. | 14. 36" Extension Tubing Set . . . . .    | 1 ea. |
| 7. Universal Disposable Speculum . . . . . | 1 ea. | 15. Disposable Cervical Dilator . . . . . | 1 ea. |
| 8. 6" Sponge Swabs . . . . .               | 2 ea. | 16. Fenestrated Drape . . . . .           | 1 ea. |

Download IFU



## PREP IN 1 STERILE FIELD

Our HSG procedure trays have the PVP solution, lubricating jelly and PVP ointment all inside the tray reducing procedure prep and set up time.

## PVC FREE


All of our HS catheters are PVC free and can be purchased separately or part of (inside) the procedure tray.

## WELCH ALLYN UNIVERSAL SPECULUM

Our HSG procedure trays contain a Welch Allyn Universal Speculum that can be used with both corded and cordless light sources.

## TECHNICAL DATA

Sterilization Method . . . . .	ETO	Speculum Type . . . . .	Welch Allyn Universal Medium
Catheter available inside procedure tray upon request		Box Qty. . . . .	10 ea.

 Not made with natural rubber latex. PVC Free & non-DEHP

## HSG PROCEDURE TRAY ORDERING INFORMATION

TMI1157 . . . . .	HSG Procedure Tray without HS Catheter
TMI1160 . . . . .	HSG Procedure Tray with 5Fr HS Catheter
TMI1161 . . . . .	HSG Procedure Tray with 7Fr HS Catheter
TMI1198 . . . . .	HSG Procedure Tray with 5Fr Shapeable Catheter
TMI1199 . . . . .	HSG Procedure Tray with 7Fr Shapeable HS Catheter



Thomas Medical is an industry leading US manufacturer of single-use women's health products.  
www.thomasmedical.com | (800) 556-0349 | sales@thomasmedical.com  
6102 Victory Way, Indianapolis, IN 46278 USA

Made in the USA  
ISO 13485 and FDA Registered

