

Simple and effective way  
to allow visibility of the  
uterine cavity



**FEATURES**

- Placement sheath decreases the need for a tenaculum
- The distal end of the catheter is balloonless allowing for a more complete view of the uterus during procedure
- An ink band provides a reference mark to indicate that proper catheter positioning has been received
- The round, closed catheter tip eases the introduction of the catheter into the uterus
- Sliding sheath assist in catheter placement

 Not made with natural rubber latex. PVC Free & non-DEHP

**TECHNICAL DATA**

French Size . . . . .	5Fr	Maximum Balloon Volume . . . . .	5Fr, 1.5ml/7Fr, 3ml
Standard Length . . . . .	28cm	Box Qty. . . . .	10 ea.
Sterilization Method . . . . .	ETO		

**ORDERING INFORMATION**

TMI1170 . . . . . Sono-Inject Catheter, 5Fr


Download IFU





Simple and effective way  
to allow visibility of the  
uterine cavity

### FEATURES

- Specially designed tip helps prevent leakage
  - Rigid catheter sheath allows for simple placement and single operator use
  - Long tip at distal end improves seating in cervical canal and has a wide tip opening for increased flow
  - Balloonless catheter allows for a more complete view of the uterine cavity during procedure
-  Not made with natural rubber latex. PVC Free & non-DEHP

### TECHNICAL DATA

French Size . . . . .	9 Fr	Sterilization Method . . . . .	ETO
Length . . . . .	30 cm	Box Qty. . . . .	10 ea.

### ORDERING INFORMATION

TMI1162-HF . . . . . eZInject SIS Catheter with High-Flow



Catheter Tip Length



Catheter Tip Opening