

THOMAS
MEDICAL

FETAL MONITORS, DOPPLERS
& ACCESSORIES



What is

THE SIGNIFICANCE OF ACCURACY?

Why buy a Toitu?

reliability
experience
advanced technology

reliability

Toitu products have proven to last many years and rarely need service.

experience

Toitu has over 30 years of experience in manufacturing high quality fetal monitors and dopplers.

advanced technology

The MT-516 has the ability detect and differentiate fetal movement of twin fetus simultaneously from mid-term to birth.



MT-516

Fetal Monitors



FD-380



FD-390



FD-491



FD-492

Fetal Dopplers

f Twin Transducer Set



Dimensions: 9.5"W x 10"H x 8"D
Weight: 9.9 lbs.

MT-516

Fetal Monitor


- a. AA2026AA Fetal Monitor "NST"
- b. AA2032AA Fetal Monitor Cart
- c. GA0109AA Y-Shaped Doppler & External UC Transducer (TR-686-05S)
- d. GA0009BA Remote Marker Switch SW-524-01 Single
- e. GA0117AA Transducer Belts, 2 ea. per set (1 ea. for doppler & 1 ea. for UC)
- f. GA1766AA Twin Transducer Set
Doppler Transducer
Personal Belt
Remote Marker Switch SW-524-02 Twin

Fetal Monitor Recording Paper
0030-022-10 Recording Paper for MT-516

The Significance of Fetal Movement

Toitu MT-516 Fetal Monitor is the only monitor that combines fetal movement with fetal heart rate and uterine contractions.

Fetal movement (FM) is shown to the right in FM waveform as well as FM indicators.



- FM Indicator (Twin 1)
- FM Indicator (Twin 2)
- FM Waveform (Twin 1)
- FM Waveform (Twin 2)

Fetal Monitor Recording Strip

FD-380

Handheld Fetal Doppler



Dimensions: 3.3"W x 6"H x 1.2"D
Weight: .5 lbs.



Dimensions: 8.5"W x 8"H x 4.5"D
Weight: 5 lbs.

FD-390

Fetal Doppler



a. BA2039AA

b. BA2042AA


c. GA0027AA

FD-380 Handheld Fetal Doppler

FD-390 Fetal Doppler

Binaural Earphone

**Compatible with FD-380 & FD-390*



FD-491

Fetal Doppler



Dimensions: 8"W x 7.5"H x 5"D
Weight: 3.5 lbs. (Includes
transducer probe & battery pack)

a. BA2039AA

b. BA2048AA

FD-491 Fetal Doppler

FD-492 Wireless Fetal Doppler

FD-490
SERIES

NEW
ARRIVAL
TO THE
TOITU LINE

FD-492

Wireless Fetal Doppler



Dimensions: 8"W x 7.5"H x 5"D
Weight: 3.5 lbs. (including
transducer probe & battery pack)



MT-516



FD-380



FD-390



FD-491



FD-492

THE SIGNIFICANCE OF ACCURACY

Be *confident* in what you are hearing

A note from the manufacturer

We would like to express our thanks to your loyal patronage.

Toitu is a medical equipment manufacturer who specializes in OB/GYN products and has started our business with the sales of Colposcope in 1958.

Ever since its foundation, we have been working on developing new products as Suction Unit, Doppler Fetus detector, Fetal Monitor and Central Station used in OB/GYN and continuing our efforts to provide attractive products which meet the needs of medical practice by advancing collaborative research with doctors and their staff. Also we have been providing perinatal management system by working hand-in-hand with medical practice since the medical information has been getting more and more important in current medical practice.

OB/GYN Department targets at treating various maladies such as perinatal, female infertility, tumors and menopause and has a broad range of areas supporting women's entire life.

We will keep on growing as "OB/GYN professionals" to make a proper suggestion on equipment and system needed for wide variety of medical practice and do a proper repair and maintenance by utilizing our longstanding relationship of trust with our OB/GYN customers and know-how we have acquired for more than half a century.

Sincerely yours,



Kazuhiro Ikeda, President

Thomas Medical, Inc. is the
proud, **exclusive distributor**
for Toitu in the USA



reliability

experience

advanced technology

THOMAS
MEDICAL

THE WORLD'S SOURCE
FOR OB/GYN PRODUCTS

a CRI[®] Company

5610 W. 82nd Street Indianapolis, IN 46278 USA
P: 800.556.0349 F: 888.411.3754
sales@thomasmedical.com

www.thomasmedical.com

