

Simple and effective way
to allow visibility of the
uterine cavity

FEATURES

- Placement sheath decreases the need for a tenaculum
- The distal end of the catheter is balloonless allowing for a more complete view of the uterus during procedure
- An ink band provides a reference mark to indicate that proper catheter positioning has been received
- The atraumatic catheter tip eases the introduction of the catheter into the uterus

TECHNICAL DATA


French Size	5Fr	Maximum Balloon Volume	5Fr, 1.5ml/7Fr, 3ml
Standard Length	31 cm	Box Qty.	10 ea.
Sterilization Method	EO		

ORDERING INFORMATION

TMI1170 Sono-Inject Catheter, 5Fr

Download IFU





Simple and effective way
to allow visibility of the
uterine cavity

FEATURES

- Specially designed tip helps prevent leakage
- Rigid catheter sheath allows for simple placement and single operator use
- Long tip at distal end improves seating in cervical canal and has a wide tip opening for increased flow
- Balloonless catheter allows for a more complete view of the uterine cavity during procedure

TECHNICAL DATA

Tip O.D.	2.0mm - 3.7mm	Sterilization Method	EO
Length	30 cm	Box Qty.	10 ea.

ORDERING INFORMATION

TMI1162-HF eZInject SIS Catheter with High-Flow



Catheter Tip Length



Catheter Tip Opening