

Enhanced features to support visualization of the uterus during pelvic laparoscopic procedures.

TMI6003

K-UMI Uterine Manipulator/Injector

IDEAL FOR

Laparoscopic Procedures, Laparoscopic Tubal Occlusion, Diagnostic Laparoscopy, Salpingoplasty


BENEFITS

Reduced catheter shaft diameter (4.5mm) allows for easy insertion

Unique pear-shaped balloon conforms to anatomy and allows for a secure seal of the cervical opening to prevent injection leaks

Flexible extension prevents cracking of luer fitting

Atraumatic beveled catheter tip reduces trauma to the uterine wall.



TMI1151

TMI6003

TMI1171

UNIQUE HANDLE CHARACTERISTICS



TMI6003 K-UMI UTERINE MANIPULATOR/INJECTOR

Pinch-grip handle allows for quick and effective repositioning during procedures



TMI1151 UTERINE MANIPULATOR/INJECTOR

Snap-on handle allows to easily place the cervical stop against the cervical wall



TMI1171 EZ GLIDE UTERINE MANIPULATOR/INJECTOR

Sliding handle adjusts after insertion and locks into place to maintain a secure atraumatic hold on the uterus

Handle features larger grip for better leverage during manipulation

ORDERING INFORMATION

TMI1151 Uterine Manipulator/Injector
 TMI6003 K-UMI Manipulator/Injector
 TMI1171 EZ Glide Manipulator/Injector

COMPLEMENTARY PRODUCTS

OS-01-025 os Finders, Box of 25
 US-01-025P Printed Uterine Sound, Box of 25
 Electrodes Various Sizes and Shapes

 Not made with natural rubber latex. PVC Free & non-DEHP

TECHNICAL DATA

UMI O.D./French Size 4.5mm/13.5Fr	Sterilization Method ETO
K-UMI O.D./French Size 5.0mm/15Fr	Maximum Balloon Volume 10cc
EZ Glide O.D./French Size 5.0mm/15Fr	Length 33cm/13"
	Box Qty. 12 ea.

Download IFU



Thomas Medical is an industry leading US manufacturer of single-use women's health products.

www.thomasmedical.com | (800) 556-0349 | sales@thomasmedical.com
 6102 Victory Way Indianapolis, IN 46278 USA

