


NEW DESIGN

Every sampler is vacuum tested for reliable performance.



FEATURES


Depth markings guide insertion depth and aid in catheter placement

Atraumatic tip reduces discomfort upon insertion

Firm and flexible cannula prevents kinking

New handle design creates a better grip during retrieval of the sample

Endometrial sampler is now non-DEHP & PVC Free



Not made with natural rubber latex.

TECHNICAL DATA

Depth Mark Range	5-11cm	100% Vacuum Tested
Length	24.5cm	French Size
Sterilization Method	ETO	Box Qty.
		25 ea.

ORDERING INFORMATION

TMI1175 Endometrial Sampler

Download IFU

