


Everything is in one sterile field to assist in an efficient and effective procedure



KIT CONTENTS

- | | | | |
|--|-------|---|-------|
| 1. Clear Tray 16"x6"x2" | 1 ea. | 9. Sanitary Napkin | 1 ea. |
| 2. Clear Tray 5"x4"x2" | 1 ea. | 10. 3"x3" Gauze Pads | 5 ea. |
| 3. Povidone-Iodine Ointment | 1 ea. | 11. 18-gauge Safety Needle | 1 ea. |
| 4. Povidone-Iodine Solution | 2 ea. | 12. 8" Swabs | 3 ea. |
| 5. Lubricating Jelly | 1 ea. | 13. 20cc Syringe | 1 ea. |
| 6. 30"x30" Over-wrap (not shown) | 1 ea. | 14. 36" Extension Tubing Set | 1 ea. |
| 7. OBP Lighted Speculum | 1 ea. | 15. Disposable Cervical Dilator | 1 ea. |
| 8. 6" Sponge Swabs | 2 ea. | 16. Fenestrated Drape | 1 ea. |

Download IFU



PREP IN 1 STERILE FIELD

Our HSG procedure kits have the PVP solution, lubricating jelly and PVP ointment all inside the kit reducing procedure prep and set up time.

PVC FREE

All of our HS catheters are PVC free and can be purchased separately or part of (inside) the procedure kit.

QUALITY TESTED


All of our HS catheters are leak-tested at both the hub and balloon.

NEW! SINGLE-USE LIGHTED VAGINAL SPECULUM

Each ER-SPEC comes complete with an integrated high powered, battery-operated LED light source.

TECHNICAL DATA

Sterilization Method	ETO	Speculum Type	ER-SPEC OBP Lighted
Catheter available inside procedure kit upon request		Speculum, Medium	
		Box Qty.	10 ea.

 Not made with natural rubber latex. PVC Free & non-DEHP

HSG PROCEDURE KIT ORDERING INFORMATION

TMI1157	HSG Procedure Kit without HS Catheter
TMI1160	HSG Procedure Kit with 5Fr HS Catheter
TMI1161	HSG Procedure Kit with 7Fr HS Catheter
TMI1198	HSG Procedure Kit with 5Fr Shapeable Catheter
TMI1199	HSG Procedure Kit with 7Fr Shapeable HS Catheter



Thomas Medical is an industry leading US manufacturer of single-use women's health products.
www.thomasmedical.com | (800) 556-0349 | sales@thomasmedical.com
6102 Victory Way, Indianapolis, IN 46278 USA

Made in the USA
ISO 13485 and FDA Registered